



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **URS-3**  
 CATALOG #: URS-3-100  
 LOT #: 27467

MANUFACTURE DATE: 09/02/2010  
 EXP. DATE: 12/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
URS Paper - Protein	URSP-PRO	27227	12/2012
URS Paper - Glucose	URSP-GLU	27067	12/2012
URS Paper - pH	URSP-PH	27311	12/2012

**TOLERANCE LEVELS OF CONTROLS**

CONTROL	LOT #	EXP. DATE	PROTEIN	GLUCOSE	pH
KOVA I	137127-5	05/2011	100 to $\geq$ 300 mg/dl	500 to $\geq$ 1000 mg/dl	7.0 – 8.0
KOVA II	138183-2	07/2011	Negative (-) to 30 mg/dl	Trace to 500 mg/dl	6.5 – 8.5
KOVA III	137752-3	02/2012	Negative (-)	Negative (-)	5.0 – 6.5

**TEST VALUES/RESULTS**

CONTROLS	PROTEIN		GLUCOSE		pH		# PASS
	Result 1	Result 2	Result 1	Result 2	Result 1	Result 2	
LEVEL I	300	300	1000	1000	7.0	7.0	6 / 6
LEVEL II	Trace	Trace	250	250	7.5	7.5	6 / 6
LEVEL III	-	-	-	-	6.5	6.5	6 / 6

COMMENTS: VISUAL CHECK OK.

APPROVED  
 BY:

*[Signature]*  
 Quality Control

DATE: 09/02/2010